

# Online Research @ Cardiff

This is an Open Access document downloaded from ORCA, Cardiff University's institutional repository: <https://orca.cardiff.ac.uk/id/eprint/134591/>

This is the author's version of a work that was submitted to / accepted for publication.

Citation for final published version:

Pennance, Tom, Archer, John, Lugli, Elena, Rostron, Penny, Llanwarne, Felix, Ali, Said M., Amour, Amour Khamis, Suleiman, Khamis Rashid, Li, Sarah, Rollinson, David, Cable, Jo ORCID: <https://orcid.org/0000-0002-8510-7055>, Knopp, Stefanie, Allan, Fiona, Ame, Shaali M. and Webster, Bonnie Lee 2020. Development of a molecular snail xenomonitoring assay to detect *Schistosoma haematobium* and *Schistosoma bovis* infections in their *Bulinus* snail hosts. *Molecules* 25 (17) , 4011. 10.3390/molecules25174011 file

Publishers page: <http://dx.doi.org/10.3390/molecules25174011>  
<<http://dx.doi.org/10.3390/molecules25174011>>

Please note:

Changes made as a result of publishing processes such as copy-editing, formatting and page numbers may not be reflected in this version. For the definitive version of this publication, please refer to the published source. You are advised to consult the publisher's version if you wish to cite this paper.

This version is being made available in accordance with publisher policies.

See

<http://orca.cf.ac.uk/policies.html> for usage policies. Copyright and moral rights for publications made available in ORCA are retained by the copyright holders.



Article

# Development of a molecular snail xenomonitoring assay to detect *Schistosoma haematobium* and *S. bovis* infections in their *Bulinus* snail hosts

Tom Pennance <sup>1,2,3,\*</sup>, John Archer <sup>1,3,\*</sup>, Elena Lugli <sup>1</sup>, Penny Rostron <sup>1</sup>, Felix Llanwarne <sup>1,4</sup>, Said M. Ali <sup>5</sup>, Amour Khamis Amour <sup>5</sup>, Khamis Rashid Suleiman <sup>5</sup>, Sarah Li <sup>6</sup>, David Rollinson <sup>1,3</sup>, Jo Cable <sup>2</sup>, Stefanie Knopp <sup>7,8</sup>, Fiona Allan <sup>1,3</sup>, Shaali M. Ame <sup>5</sup> and Bonnie L. Webster <sup>1,3</sup>

<sup>1</sup> Wolfson Wellcome Biomedical Laboratories, Department of Zoology, Natural History Museum, Cromwell Road, London SW7 5BD, UK; e.lugli@nhm.ac.uk (E.L.); p.rostron@nhm.ac.uk (P.R.); felix.llanwarne1@student.lshnm.ac.uk (F.L.); d.rollinson@nhm.ac.uk (D.R.); f.allan@nhm.ac.uk (F.A.); b.webster@nhm.ac.uk (B.L.W.)

<sup>2</sup> School of Biosciences, Cardiff University, Cardiff, CF10 3AX, UK; cablej@cardiff.ac.uk

<sup>3</sup> London Centre for Neglected Tropical Disease Research (LCNTDR), London, UK.

<sup>4</sup> Faculty of Infectious and Tropical Diseases, London School of Hygiene and Tropical Medicine, Keppel Street, London, WC1E 7HT, UK

<sup>5</sup> Public Health Laboratory – Ivo de Carneri, P.O. Box 122 Chake-Chake, Pemba, United Republic of Tanzania; saidmali2003@yahoo.com (S.M.A.); amourkhamis2003@yahoo.com (A.K.A.); khasule66@yahoo.com (K.R.S.); shaaliame@yahoo.com (S.M.A.)

<sup>6</sup> Schistosomiasis Resource Centre, Biomedical Research Institute, 9410 Key West, Rockville, MD 20850, United States of America; sli@afbr-bri.org

<sup>7</sup> Swiss Tropical and Public Health Institute, Socinstrasse 57, 4002 Basel, Switzerland; s.knopp@swisstph.ch

<sup>8</sup> University of Basel, Petersplatz 1, 4001 Basel, Switzerland

\* Correspondence: t.pennance@nhm.ac.uk (T.P.); j.archer@nhm.ac.uk (J.A.); Tel.: 020-7942-6884 (T.P.)

Received: 31 July 2020; Accepted: date; Published: date

**Abstract:** Schistosomiasis, a neglected tropical disease of medical and veterinary importance, transmitted through specific freshwater snail intermediate hosts, is targeted for elimination in several endemic regions in sub-Saharan Africa. Multi-disciplinary methods are required for both human and environmental diagnostics to certify schistosomiasis elimination when eventually reached. Molecular xenomonitoring protocols, a DNA based detection method for screening disease vectors, have been developed and trialed for diseases transmitted by hematophagous insects, such as lymphatic filariasis and trypanosomiasis, yet none have been extensively trialed or proven reliable for schistosomiasis. Here, previously published universal and *Schistosoma* specific internal transcribed spacer (ITS) rDNA primers have been adapted into a triplex PCR primer assay that allows for simple, robust and rapid detection of *Schistosoma haematobium* and *S. bovis* in *Bulinus* snails. We show this two-step protocol can sensitively detect DNA of a single larval schistosome from experimentally infected snails and demonstrate its functionality for detecting *S. haematobium* infections in wild caught snails from Zanzibar. Such surveillance tools are a necessity for succeeding in and certifying the 2030 control and elimination goals set by the World Health Organization.

**Keywords:** bovine; control; elimination; schistosomiasis; urogenital; surveillance; disease; parasite

## 1. Introduction

Schistosomiasis is a parasitic disease infecting an estimated 229 million people worldwide caused by parasitic worms of the genus *Schistosoma*, infection leading to severe morbidity and

mortality due to the associated complications of worm presence [1]. *Schistosoma* spp. in Africa are transmitted through specific freshwater snail intermediate hosts of the *Bulinus* and *Biomphalaria* genus [2]. Infections occur when humans or animals come into contact with freshwater containing infectious larval stages (cercariae) shed from the infected snails. Human schistosomiasis in Africa, where at least ~90% of those requiring treatment live [3], consists of two forms of the disease, urogenital and intestinal schistosomiasis, caused predominantly by *Schistosoma haematobium* and *S. mansoni* respectively [1]. Bovine, ovine and caprine schistosomiasis is also of significant veterinary and economic importance across sub-Saharan Africa [4,5] and is caused by infection of cattle, sheep and goats with species closely related to *S. haematobium* (termed *S. haematobium* group species), primarily *S. bovis*, *S. curassoni* and *S. mattheei*. Sympatric geographical distribution of schistosome species and overlapping intermediate snail host use of certain schistosome species and geographical strains, complicates disease transmission surveillance in (co)endemic zones [2,6,7].

The World Health Organization (WHO) aims for the elimination of human schistosomiasis as a public health problem, defined as >1% of the population with heavy intensity infections ( $\geq 50$  schistosome eggs per 10 ml of urine or  $\geq 400$  schistosome eggs per gram of faeces [8]), in all endemic countries by 2030 [9]. Despite great advances in schistosomiasis control mainly via preventative chemotherapy (praziquantel), the lack of protection against rapid re-infection together with the prolific asexual replication of schistosomes within their intermediate snail host presents substantial hurdles to achieving the targeted elimination of schistosomiasis. Very quickly snails can become infected by eggs emanating from untreated humans leading to rapid resurgence of transmission [10]. Therefore, adaptive treatment strategies that take into account the transmission dynamics of *Schistosoma* spp. with their snail hosts are required to control and eliminate the disease [11].

To better understand the local transmission dynamics of different *Schistosoma* species, allowing both human and bovine schistosomiasis to be monitored, a need exists for methodologies that detect schistosome infections in the intermediate host snails. These tools for assessing *Schistosoma* transmission could eventually be used during elimination programs to identify focal areas of persisting transmission or certify elimination and / or transmission interruption [12–14]. Defining ongoing transmission in snail populations through traditional methods of observing cercariae shed from snails is particularly challenging in an elimination setting such as the Zanzibar Archipelago where few snails (0.5 – 2.3%) are observed shedding cercariae [6,15]. Furthermore, snails with non-patent (including pre-patent) infections are missed using these approaches. Where cercarial larval schistosomes are observed, these cannot be easily identified using morphological characteristics to species level (although relative position of sensory receptors is of some value [16,17]).

Molecular xenomonitoring is a DNA-based method that has been developed to monitor the transmission of several vector borne diseases including trypanosomiasis [18,19], filariasis and malaria [20], helminthiasis [21] and fascioliasis [22], including to some extent schistosomiasis [23–29]. Screening snails provides evidence on the extent of environmental contamination (i.e. schistosome miracidia penetrating snails) as well as environmental infection risk (i.e. schistosome sporocysts and cercariae developing inside the (pre-patent) snails, eventually emerging from the snail (patent)). Most of the available snail-schistosome xenomonitoring assays do not include internal controls [23,28,30], an important feature in any diagnostic tool that helps prevent false negative results [27]. Many assays just assume that a negative result means non-infection.


In the current study, we adapted available universal [31] and *Schistosoma* specific [27] ITS rDNA primers to design a three primer multiplex assay and tested this as a simple, robust and rapid xenomonitoring PCR assay to enable the large-scale screening of *Bulinus* snails for *Schistosoma* infections (*S. haematobium* and *S. bovis*). We use a conventional PCR based approach focused on simplicity, ease of data interpretation, sensitivity and specificity with a primary aim to provide a xenomonitoring tool for monitoring *S. haematobium* transmission in endemic settings.

## 2. Results

### 2.1. In silico and in vitro Primer Evaluation

*B. globosus* and *B. nasutus* rDNA sequence data showed conserved primer binding sites for the universal primers ETTS2 and ETTS1 [31] at the 3' end of the 18S and 5' end of the 28S, flanking regions of the ITS, respectively. ETTS1 gave a 100% match and the ETTS2 primer showed just a single base pair mismatch. The resulting snail amplicon size predicted from these alignments was between 1232–1263 bp and served as an internal snail control during PCR amplification.

Alignments of the *Schistosoma* specific ITS primers (ITS2\_Schisto\_F and ITS2\_Schisto\_R [27]) showed 100% and 90% (2 mismatches) homology to *S. haematobium* and *S. bovis* respectively, with no cross reactivity to the *Bulinus* reference rDNA data. When paired with their opposing universal primers (ITS2\_Schisto\_F + ETTS1 or ITS2\_Schisto\_R + ETTS2) amplicon sizes of 538 and 835 bp were predicted respectively for *Schistosoma*. With the addition of the other universal primer to each combination (ETTS2 and ETTS1 respectively), the three-primer Multiplex ITS Xenomonitoring (MIX) reactions were predicted to be able to produce distinct amplicon profiles for non-infected snails (a single snail amplicon) and snails infected with *Schistosoma* spp. (three band profile). This was confirmed by in vitro testing of the primer combinations (Figure 1).



**Figure 1.** Singleplex (A; ETTS2 + ETTS1) and multiplex (B; multiplex ETTS2 + ETTS1 + ITS2\_Schisto\_F, C; ETTS2 + ETTS1 + ITS2\_Schisto\_R) PCRs on laboratory bred *Bulinus wrighti* (*B.w.*) and *Schistosoma haematobium* (*S.h.*) gDNA separately (1; *B.w.*, 2; *S.h.*) and combined (3; *B.w.* + *S.h.*). When *B.w.* and *S.h.* DNA is combined (A3, B3, C3), two amplicons are produced by the ETTS1 + ETTS2 primers, a larger snail amplicon (Sn) (~1200 bp) and a smaller *Schistosoma* amplicon (T) (~1000), with the additional *Schistosoma* specific primers producing either a 538 bp (B3; ITS2\_Schisto\_F) or 835 bp (C3; ITS2\_Schisto\_R) amplicon (S). A larger amplicon (A) (~1400–1600 bp) was also observed to be amplified in some reactions, this is thought to be a PCR artefact or additional primer targets in the *Schistosoma* gDNA. L1 = HyperLadder I (Bioline, London, UK). -ve = negative no template control.


To maximise amplification efficiency/sensitivity and to provide good amplicon size differentiation, the multiplex PCR incorporating the internal ITS2\_Schisto\_F (Figure 1B) was selected for further development and testing. This primer combination was also selected as it targets the ITS2 region for *Schistosoma* containing four species specific SNP's enabling species identification (Table 1).

The MIX assay proved robust at varying annealing temperatures (55°C, 60°C Figure 2A, 58°C Figure 2B) and with 58°C proving to be the most efficient, maximising specificity without decreasing sensitivity. Each of the three amplicons were extracted from the gel and sequenced, confirming the band identity and specificity of the primers to their target gDNA amplicon. These three bands are described as the snail (Sn) (1232 – 1263 bp), trematode (T) (~ 1000 bp) and *Schistosoma* (S) (538 bp) bands going forward. The secondary *Schistosoma* ITS xenomonitoring (SIX) PCR, solely targeting the *Schistosoma* amplicon proved robust, enabling single amplicon generation and sequencing (Figure 3). This provides a two-step PCR methodology with the MIX PCR for the initial high-throughput screening of the samples and the secondary SIX PCR to target specific samples for further infection clarification.




**Table 1.** *Schistosoma* species-specific SNP positions (including base position) in the ITS2 region.

<i>Schistosoma</i> species	ITS 2 <i>Schistosoma</i> species specific SNP positions (bp)			
	SNP1 (90)	SNP2 (145)	SNP3 (195)	SNP4 (265)
<i>S. haematobium</i>	<i>S. h</i> (G)	<i>S. h</i> (C)	<i>S. h</i> (G)	<i>S. h</i> (C)
<i>S. bovis</i>	<i>S. b</i> (A)	<i>S. b</i> (T)	<i>S. b</i> (A)	<i>S. b</i> (T)




**Figure 2.** Multiplex ITS xenomonitoring assay trial at 55°C (A) and 60°C (B). Includes gDNA of *Bulinus wrighti* of both BioSprint (Lane 1 and 5) and DNeasy extractions (Lane 2 and 6) and gDNA of *Schistosoma haematobium* (Lane 3 and 5) and *S. bovis* (Lane 4 and 6). Combinations of *B. wrighti* and *S. haematobium* (Lane 5) or *S. bovis* (Lane 6) gDNA shown. Sn = snail amplicon, T = trematode amplicon, S = *Schistosoma* amplicon and A = non-specific amplicon or artefact. L1 = HyperLadder I. L2 = HyperLadder IV (Bioline, London, UK). -ve = negative no template control.



**Figure 3.** Gel showing the secondary singleplex ITS xenomonitoring (SIX) PCR for: 1) *Schistosoma haematobium* gDNA; 2) *S. bovis* gDNA; 3) *S. haematobium* + *B. wrighti* gDNA; 4) *S. bovis* + *B. wrighti* gDNA. -ve = non-template negative control. L1 = HyperLadder I (Bioline, London, UK).

## 2.2. Analytical Sensitivity

The assay proved highly sensitive with a LoD of 0.02 ng and 0.002 ng of gDNA for *S. bovis* and *S. haematobium*, respectively (Figure 4). Sensitivity appeared higher for *S. haematobium* (Figure 4), but in both cases the assay's sensitivity is above that necessary to detect gDNA from a single miracidium, which ranges from 1.6–3.65 ng/μl [32]. At lower *Schistosoma* DNA concentrations the 1005 bp trematode band (T) lost sensitivity compared with the smaller *Schistosoma* specific band.




**Figure 4.** Sensitivity tests of ITS1-2-F PCR performed with serial dilutions of *Schistosoma haematobium* and *S. bovis* gDNA in the presence of *Bulinus wrighti* gDNA. L1 = HyperLadder I. L2 = HyperLadder IV (Bioline, London, UK).

### 2.3. Experimental snail infections

For the non-patent infections of *B. truncatus* with *S. haematobium*, preserved 24 h after exposure, 61.1% (11 out of 18) *B. truncatus* were observed to be infected presenting the *Schistosoma* specific ITS2 band (Figure 5). Infections were detected in snails exposed to 1, 2 and 7 miracidia. Two of the five (40%) *B. truncatus* exposed to one or two miracidia and left for 11 weeks, did not reach patency, but were also confirmed to be infected by the presence of *S. haematobium* amplicon (Figure 5: Lanes 20 and 21). The secondary SIX PCR was performed on all 13 non-patent infected snails and the single amplicons were sequenced and confirmed as *S. haematobium*. Out of all the snails infected that survived until the end of the experiment (11 weeks), 15% (nine out of 62) reached patency of which two had been infected with two miracidia and seven with seven miracidia. One of these samples, infected with two miracidia, was analysed using the MIX PCR giving the expected triple banding pattern (snail, trematode and *Schistosoma*) (Figure 5: Lane 24). All three amplicons from this sample (Figure 5: Lane 24) were gel extracted and sequenced confirming their identification. Interestingly, in all the non-patent infections the large trematode amplicons (ETTS2-ETTS2) did not amplify (Figure 5), due to the low level of *Schistosoma* DNA present in the snails that did not reach patency.



to amplify but the remaining 31 produced the *Schistosoma* amplicon that all sequenced as *S. haematobium*. One sample also gave the trematode band without the *Schistosoma* band indicating a non-*Schistosoma* trematode infection.



**Figure 6.** Gel images for the multiplex ITS xenomonitoring (MIX) PCR amplicons for 94 non-patent *Bulinus globosus* collected from Wambaa, Pemba, United Republic of Tanzania. Text under each amplicon denotes the outcome of the *Schistosoma* sp. targeted sequencing where relevant (i.e. presence of 538 bp amplicon), which resulted in either *S. haematobium* (S.h.) or sequencing failure (F). Presence of a trematode band without the presence of the *Schistosoma* band indicates a non-*Schistosoma* trematode infection (Tr.). Other non-specific bands, in this case larger bands (NA) were also observed in these snail populations, which did not amplify with the secondary SIX PCR. x = sample failure with no control amplicon. Arrows highlight the presence of the ~1000 bp trematode band when present (n = 8). *B. globosus* with a patent *S. haematobium* infection (Cham10.1 see [6]) was run as a positive control (+ ve) and also represents the amplicon profile obtained for the seven patent *B. globosus* snail (five and two with *S. bovis* and *S. haematobium* infections, respectively (see section 2.4). -ve = the non-template negative control. L1 – HyperLadder IV (Bioline, London, UK).

## 2.6. *Schistosoma* spp. *cox1* RD-PCR

Despite trying different annealing temperatures and gDNA template amounts used, the *cox1* RD-PCR developed by Webster *et al.* (2010) [33] tested on the patent *S. haematobium* and *S. bovis* infected *B. globosus*, only generated the species-specific amplicon for *S. bovis* infected snails. PCRs for the *S. haematobium* infected snails repeatedly failed to produce a clear amplicon. The *cox1* amplicons produced for the *S. bovis* infected snails were sequenced and the data matched that obtained from the cercariae collected and analysed from the snails (see [6]).

## 3. Discussion

Pre-patent and non-patent snail screening methods for schistosomes, such as molecular xenomonitoring, offer a higher sensitivity over traditional snail shedding methods that can only detect patent infections (observation of schistosome cercariae). Molecular xenomonitoring better helps to assess the impact of schistosomiasis control interventions in local communities, particularly where local elimination is being achieved and certification of the absence of transmission is suspected at specific foci. However, the diversity of the *Schistosoma* species circulating in co-endemic areas means that species specific methodologies are needed to prevent false positive data related to non-target species.



Here, we describe the development and application of a molecular xenomonitoring pipeline for the detection and differentiation of *S. haematobium* and *S. bovis* patent and non-patent infections in *Bulinus* freshwater snails, using three previously developed primers [27,31]. The MIX assay screens for *Schistosoma* and other trematode species, whilst also incorporating an internal control, in this case gastropod DNA, an important feature for any molecular diagnostic assay. The MIX PCR generates clearly identifiable amplicons, of different sizes, for each target (snail, trematode, *Schistosoma*) which are visible by simple agarose gel electrophoresis. However, the trematode target lacks sensitivity at low DNA concentrations, probably due to its large size and PCR biases for small amplicons at reduced gDNA concentrations. Interestingly, a PCR artefact (~1400-1600 bp) was also observed when using the MIX assay in the presence of *Schistosoma* DNA, suggesting that the primers may have a secondary binding site. However, this artefact is clearly identifiable from the main target amplicons and does not mislead interpretation of the results.

### 3.1. Sensitivity of MIX PCR assay

Our *in silico* and *in vitro* testing of the MIX assay showed that the presence of *S. haematobium* and *S. bovis* DNA can be routinely detected at low concentrations, and also was able to identify non-patent *Schistosoma* infections in snails where the level of DNA varies depending on the development of the infection. The LoD for *Schistosoma* DNA was  $\leq 0.02$  ng/ $\mu$ l, which is 80-fold lower than the minimum amount of gDNA usually observed from a single miracidia [32]. This was also demonstrated by the assay's ability to detect pre-patent snail infections 24 h after exposure to a single miracidium. This provided sufficient sensitivity for the LoD of detection needed to detect any stage of snail infection, from initial miracidial penetration of a single miracidium to full patency, in natural settings. The fact that not all the snails tested from the experimental snail infections gave positive results is corroborative with observations that, even in experimental systems, many snails avoid penetration or destroy the miracidia rapidly upon invasion. The MIX and SIX methodology also proved robust when used to screen 'wild caught' snails from Pemba, with uninfected, pre-patent *S. haematobium* infected snails, and non-*Schistosoma* trematode infections clearly identified.

### 3.2. Benefits of an updated molecular xenomonitoring protocol for schistosomiasis surveillance

The molecular xenomonitoring protocol requires few consumables and no cold chain, and results can be interpreted using basic molecular laboratory equipment (thermocycler and gel electrophoresis) making the molecular assay accessible in lower resource settings, such as schistosomiasis endemic regions. The molecular xenomonitoring approach described here therefore provides a useful tool for monitoring schistosomiasis transmission, as has been outlined as a necessary method for leading toward the WHO 2030 goals for schistosomiasis control and elimination [9].

Molecular xenomonitoring surveillance techniques are often associated with diseases transmitted by hematophagous insects, such as lymphatic filariasis in mosquito vectors [20,34–36] and trypanosomes in tsetse flies [18,19]. However, several assays have been developed for detecting trematode species in freshwater snails, including *Fasciola* spp. [22,37–44], other wildlife trematode species [45] and medically important schistosome species; *S. japonicum* [46,47], *S. mansoni* [24,27,28,48–54] and *S. haematobium* [23,26,27,30,52,54,55]. The first developed assay for the molecular detection of *S. haematobium* DNA in *Bulinus* employed the highly repetitive *Dra1* target and this has been the marker of choice for studies investigating *S. haematobium* infections in snails due to its high sensitivity [55]. However, the specificity of the *Dra1* and interpretation of results can be problematic due to the frequent false positive and negative results, lack of internal control, and difficulties in interpreting the amplicon patterns. Furthermore this marker does not allow for species identification. Kane *et al.* (2013) [54] employed the use of another repetitive marker, intergenic spacer (IGS), for the detection of snail infections and a post amplification restriction digest allowed for downstream species identification of *S. haematobium* and *S. bovis*. However the method lacks internal controls. In addition, many of these assays use quantitative-PCR (qPCR), rather than conventional PCR/gel electrophoresis. Although able to quantify levels of DNA within a sample, qPCR is more arduous to

carry out and lesser suited for use in endemic settings. However, recent technological advances in sample preparation and DNA extraction methods have demonstrated robust field setting methodologies to conduct qPCR analysis capable of detecting avian trematodes and host species in Canadian lakes [56–58], which could potentially be modified to suit the detection of human and bovine schistosomes in sub-Saharan Africa, although cost and throughput would need to be considered.

A recent assay designed by Schols *et al.* (2019) [27] is a six primer multiplex PCR, that incorporates an internal snail control and offers a xenomonitoring tool for *S. haematobium* group species that are transmitted by *Bulinus* snail hosts. Our study simplifies the multiplex process, reducing the primer numbers and mitigating against PCR competition and biases that may occur with multiple primer combinations. It also allows for greater amplicon size differentiation (as amplicon sizes can be more easily distinguished based on size) making results more interpretable. The ITS rDNA is a favourable target within the repeat ribosomal operon of *Bulinus* and *Schistosoma* spp., easily detected within small quantities of DNA due to the high copy number of rRNA clusters within eukaryote genomes. The other key feature of the target relates to specificity. The ITS regions of *Schistosoma* and *Bulinus* spp. can be routinely amplified using conventional PCR thanks to its small size ~1000 bp) and highly conserved flanking regions (5'18S and 3'28S) enabling the use of universal primers (ETTS1+2) for multiple species [31]. However inter species heterogeneity, and to a lesser extent intra species heterogeneity (Pennance *et al.*, unpublished observations), of the ITS regions allow for differentiation between species, such as those of the *S. haematobium* group [7,33]. The internal *Schistosoma* specific primer is situated in a conserved ITS region within the *Schistosoma* genus, with 100% conservation between several African species suggesting that it could be utilised for several *Schistosoma*-snail transmission systems.

### 3.3. Limitations of molecular xenomonitoring approaches for schistosomiasis surveillance

From our study, we identified two limiting factors for the practical use of the methodology. First, the laborious nature of testing each individual snail adds time and cost. Further sensitivity testing should be performed to support the development of pooling strategies. This would help to determine whether infections are still detected when the *Schistosoma* DNA is diluted in the presence of much higher concentrations of snail DNA, which may inhibit the reaction. Pooling strategies have been successful for arthropod xenomonitoring protocols [18] and would allow for higher throughput of samples required for screening large snail populations, such as those encountered for schistosomiasis.

Second, a limitation does come with the need for the secondary screening (SIX PCR) of the *Schistosoma* amplicon, via sequencing, to confirm species. Despite best efforts, rapid species diagnostics, such as the rapid diagnostic *cox1* RD-PCR developed by Webster *et al.* (2010) [33] to determine adult worm and larval stage species identity, was not robust when snail DNA was present, particularly for *S. haematobium* infections. The *cox1* RD-PCR was suggested as a secondary screening method by Schols *et al.* (2019) [27] but it was only theoretically examined as part of that study. Clearly further 'wet lab' testing on infected snails is needed. In regions where *Schistosoma* hybridisation occurs, mitochondrial DNA analysis would be necessary, since both nuclear and mitochondrial DNA is required for hybrid identification [33]. Unfortunately, as with most diagnostics there is a balance between sensitivity and specificity, with sensitivity increasing and specificity decreasing, usually due to the nature of the biomedical targets. Here, rapid screening with high sensitivity was a priority due to the extremely low levels of infections in our study sites, with secondary species-specific screening only required on a small subset of samples that were identified as infected. Moreover, Zanzibar was previously thought to be a zone of *S. haematobium* transmission only, although with the recent occurrence of *S. bovis* transmission being observed [6], the additional species specific screening is warranted. However, the need for the secondary screening step for *Schistosoma* species identification does need further exploration such as trialing more direct methods that mitigate DNA sequencing, for example, amplicon enzyme restriction digestion demonstrated in Kane *et al.* (2013) [54]. However, it is also important to gather detailed information, as is obtained through DNA sequencing and

analysis, on species complexities and diversity within target endemic zones to optimise focal surveillance strategies [6,7]. It is likely that xenomonitoring methods may need to be adapted to specific endemic zones due to geographical genetic difference of the target organisms and potential unidentified species.

#### 4. Materials and Methods

##### 4.1. Primer selection and in silico evaluation


The universal primer pair, ETTS2 and ETTS1 (Table 2 and Figure 7), was selected for the development of the internal control for the assay. They anneal to conserved flanking regions either side of the ITS(1+2) rDNA region of *Schistosoma* spp., amplifying the full ITS rDNA regions resulting in an amplicon of ~ 1005 bp [6,7,31,33]. These primers have also been demonstrated to amplify the full ITS rDNA region of other organisms including intermediate gastropod hosts. Primer cross reactivity with the target *Bulinus* snail hosts was further confirmed through alignments of the ETTS2 and ETTS1 primers with *B. globosus* and *B. nasutus* rDNA regions, available from ongoing projects (Briscoe *et al.* unpublished data, Pennance *et al.* unpublished data).

To develop the *Schistosoma* specific target, two *Schistosoma* specific primers (ITS2\_Schisto\_F and ITS2\_Schisto\_R) published by Schols *et al.* (2019) [27] were selected targeting the internal ITS1 and ITS2 rDNA regions of *Schistosoma* (Figure 7). These were further tested *in silico* for specificity by stringently aligning them with rDNA sequence data (Briscoe *et al.*, unpublished data; Pennance *et al.*, unpublished data) of a single *B. globosus* and *B. nasutus* from both Unguja and Pemba island (Zanzibar, United Republic of Tanzania) and those previously published for *Schistosoma* spp. [59,60].

All alignments were performed using Sequencher v5.4.6 (Gene Codes Corporation, Michigan, USA) and primer positions were used to predict the specific amplicon sizes that would result following amplification of snail and schistosome DNA using the different primer combinations of ETTS1, ETTS2, ITS2\_Schisto\_F and ITS2\_Schisto\_R.

**Table 2.** Details of the primers selected for the development of the xenomonitoring assay. Universal (U) and specific (S) denotes whether the primers universally targets both *Schistosoma* and snail or just specifically target *Schistosoma* DNA.

Primer (direction)	Primer Sequence (5'-3')	Primer position	State	Reference
ETTS1 (Reverse)	TGCTTAAGTTCAGCGGG	28S 5' end (ITS2 3' flanking region)	U	Kane <i>et al.</i> (1994) [31]
ETTS2 (Forward)	TAACAAGGTTTCCGTAGGTGA	18S 3' region (ITS1 5' flanking region)	U	Kane <i>et al.</i> (1994) [31]
ITS2_Schisto_F (Forward)	GGAAACCAATGTATGGGATTATTG	ITS1 3' end (5.8S 5' flanking region)	S	Schols <i>et al.</i> (2019) [27]
ITS2_Schisto_R (Reverse)	ATTAAGCCACGACTCGAGCA	ITS2 (middle)	S	Schols <i>et al.</i> (2019) [27]



**Figure 7.** Primer annealing positions flanking and internal to the ITS1+2 rDNA targets. Primer positions are mapped to *Schistosoma haematobium* and *S. bovis* ITS1+2 reference data [59], and to a *Bulinus globosus* and *B. nasutus* DNA reference (Pennance et al., unpublished data). For *Schistosoma* DNA the primer combinations produce two fragments; 1) ETTS2-ETTS1 (1005 bp) and either 2) ITS2\_Schisto\_F-ETTS1 (538 bp) or 3) ITS2\_Schisto\_R-ETTS2 (835 bp). For *Bulinus* DNA the primer combinations produce one fragment ranging in size between 1232-1263 due to inter species variation. For *Schistosoma* species identification four SNPs are present at bp positions 90, 145, 195, and 265 in the ITS2 rDNA region, allowing differentiation of *S. haematobium* and *S. bovis* following ITS2 sequencing.

#### 4.2. *Bulinus* and *Schistosoma* genomic DNA extractions

Whole soft tissue from *Bulinus* samples (as detailed below) available through the Schistosomiasis Collection at the Natural History Museum (SCAN) [61] and other ongoing projects, including laboratory and field samples, infected / non-infected and patent / non-patent, were used for the assay development and validation. Genomic DNA (gDNA) from all *Bulinus* samples were extracted using a modified tissue lysis protocol [6]. Two kits were then used to extract total gDNA from the lysed snail tissue, the BioSprint 96 DNA Blood Kit (Qiagen, Manchester, UK) for high-throughput multiple sample processing, and the DNeasy Blood & Tissue Kit (Qiagen, Manchester, UK) for single sample processing. Protocols were carried out according to the manufacture instructions.

Positive control *Schistosoma* gDNA was obtained from adult worms, *S. haematobium* (single female worm from Zanzibar) and *S. bovis* (single male worm from Senegal), available from SCAN. DNA was extracted following the DNeasy Blood & Tissue Kit protocol according to manufacturer's instructions (Qiagen, Manchester, UK) [60].

#### 4.3. PCR conditions, amplicon visualisation and sequencing

All PCR amplifications were performed in 25 µl PCR reactions using illustra™ PuReTaq Ready-To-Go™ PCR Beads (GE Healthcare, UK) with 1 µl of each primer, in their different combinations as stated in each section, at a concentration of 10 µM. gDNA templates (*Schistosoma* and/or *Bulinus* sp.) were added at different volumes and concentrations as detailed below. The PCR cycling conditions for all multiplex and singleplex reactions were as follows: initial denaturation 5 minutes at 95°C followed by 40 cycles of 30 seconds at 95°C, 30 seconds at 58°C (unless stated otherwise), 90 seconds at 72°C and a final extension of 10 minutes at 72°C. Visualisation of all PCR products were performed by running 7.5 µl of each PCR product, mixed with 2 µl of Bioline 5x DNA Loading Buffer Blue (London, UK) and GelRed for visualisation under UV light, on a 2% agarose gel for 90 minutes at 90



V. HyperLadder I and HyperLadder IV were run alongside the PCR amplicons to assess fragment sizes. Gels were visualised using a GBOX-Chemi-XRQ gel documentation system (Syngene, Cambridge, UK).

To validate amplification specificity, selected PCR amplicons from multiplex PCRs, where multiple amplicons are present, were cut from agarose gels and sequenced following purification using the QiaQuick Gel purification kit (Qiagen, Manchester, UK) following manufacturer's instructions. For singleplex reactions, resulting in a single amplicon, PCR products were purified using the QiaQuick PCR purification kit (Qiagen, Manchester, UK) following manufacturer's instructions. Amplicons were Sanger sequenced in both directions using dilutions of the PCR primers. Sequence data was manually edited using Sequencher v5.4.6 (Gene Codes Corporation, Michigan, USA) and amplicon identification was confirmed by comparison to *Schistosoma* reference data [59] and by BLAST analysis (BLAST: Basic Local Alignment Search Tool, NCBI).

#### 4.4. In vitro primer testing and assay validation

All gDNA extractions from laboratory-bred *Bulinus wrighti* (not exposed to any trematodes and therefore negative for infection) and from the *S. haematobium* and *S. bovis* adult worms were quantified using a Qubit® Fluorometer using the dsDNA Broad Range (BR) Assay Kit (Molecular Probes, Life Technologies). The gDNA extracts from the single adult *S. haematobium* and *S. bovis* worms were normalised, using nuclease free water, to 2 ng/μl (+/- 0.05 ng/μl). The gDNA extract of a *B. wrighti* snail control was recorded and kept at 31.3 ng/μl. Template gDNA (1 μl) was used in each PCR separately or combined and used to test the different primer combinations (shown in Figure 1). The primers were tested as singleplex PCRs for the internal control (ETTS2 + ETTS1) targeting both snail and *Schistosoma* gDNA and then as multiplex PCR's incorporating each of the internal *Schistosoma* specific primers (ETTS2 + ETTS1 + ITS2\_Schisto\_F or ITS2\_Schisto\_R). All test PCRs were initially performed at an annealing temperature of 55°C.

The multiplex primer combination ETTS2 + ITS2\_Schisto\_F + ETTS1 was selected and taken forward for further development and validation. This is referred to as the Multiplex ITS Xenomonitoring (MIX) PCR. The MIX PCR was further tested at annealing temperatures of 58° and 60°C to enhance assay specificity, with 58°C taken forward for further experiments. Additionally, a secondary *Schistosoma* ITS xenomonitoring (SIX) PCR, incorporating just the *Schistosoma* specific primer (ITS2\_Schisto\_F) and its universal reverse primer (ETTS1), was validated targeting just the 538 bp *Schistosoma* DNA amplicon. The SIX PCR was developed to obtain more targeted schistosome species data amplicon sequence analysis, of positive samples, following initial high-throughput screening of snail populations with the multiplex PCR, which incorporates the internal snail control.

#### 4.5. Sensitivity testing

Analytical sensitivity and limit-of-detection (LoD) of the MIX PCRs ability to detect low levels of *Schistosoma* DNA, was performed using serial dilutions of *S. haematobium* and *S. bovis* gDNA. The *S. haematobium* and *S. bovis* gDNA, normalised to 2 ng/μl (+/- 0.05 ng/μl), was diluted using nuclease free water by one in ten (0.2 ng/μl), one in one hundred (0.02 ng/μl) and one in one thousand (0.002 ng/μl). 1 μl of each *Schistosoma* gDNA dilution was used in each multiplex PCR together with 1 μl of the *B. wrighti* gDNA (31.3 ng/μl).

Sensitivity was also tested using controlled laboratory snail infections. Infections were performed by the Schistosomiasis Resource Centre (SRC) (Biomedical Research Institute, Maryland, USA [62]) using their *B. truncatus* / *S. haematobium* (Egyptian strain) model lifecycle system. Juvenile *B. truncatus* (2-3 mm, n=133) were divided into three groups, with individual snails in each group being exposed to either 1, 2 or several (~7) *S. haematobium* miracidia respectively (Table 3). Miracidia, hatched in freshwater from eggs collected from *S. haematobium* infected male LVG Syrian golden hamsters (see Ethical Statement), were added to individual *B. truncatus* snails which had been placed in fresh snail water, in individual wells of a 24 well ELISA plates. A fine tipped Pasteur pipet was used under a dissection microscope to capture and deliver either an individual miracidium or several

(~7) miracidia at a time, following the standard operating procedures (SOPs) conducted at SRC (see: <https://www.afbr-bri.org/schistosomiasis/standard-operating-procedures/>).

The snails were kept in their individual wells until no miracidia were observed swimming under a binocular microscope, assumed to have penetrated the snail (~2 h). Following 24 h after initial exposure to the miracidia, half of each infection group were preserved in 100% ethanol for molecular analysis. The remaining exposed *B. truncatus* were maintained in their separate infection groups for 11 weeks to allow the infections to mature and since this was the first opportunity to conduct sampling of infected snails. Snails were maintained according to the SRC's SOP's (see above). Snails that died were recorded and promptly removed from the group. At 11 weeks post exposure the remaining snails were individually induced to shed cercariae by exposure to fresh water and light. Once it had been determined if the snails were infected and patent they were washed, to remove any cercariae, and preserved in 100% ethanol for molecular analysis.

The MIX PCR was performed using gDNA (1 µl) extracted from six individual *B. truncatus* from each group that were preserved after 24 h, two non-patent snails from group 1 and 2, and one non-patent snail from group 3 (11 weeks post exposure), and one patent (shedding) snail from group 2 (11 weeks post exposure) (Table 3). The secondary SIX PCR was performed on selected *Schistosoma* positive samples, to amplify the 538bp *S. haematobium* specific amplicon for sequence analysis to confirm that the MIX PCR was not a false positive.

**Table 3.** Groups of *Bulinus truncatus* (*B.t.*) experimentally challenged with either 1, 2 or ~7 *S. haematobium* (*S.h.*) miracidia and preserved 24 hours (h) post exposure or checked for patent *S.h.* infections and preserved 11 weeks (wks) post exposure.

Infection Group	No. of <i>B.t.</i> exposed	No. of <i>S.h.</i> miracidia used	No. of <i>B.t.</i> preserved at 24 h	No. of <i>B.t.</i> checked for patency at 11 wks and preserved (no. shedding +ve)
1	45	1	22	22 <sup>1</sup> (0)
2	43	2	21	19 <sup>1</sup> (2)
3	45	~7	23	21 <sup>1</sup> (7)

<sup>1</sup> One *B. truncatus* died from each infection group during the 11 weeks post miracidia exposure.

#### 4.6. Specificity testing and validation on field samples

As part of a longitudinal xenomonitoring project on Pemba in relation to urogenital schistosomiasis transmission [6], 'wild caught' *B. globosus* and *B. nasutus* field isolates were available for further validation of the MIX assay. Individual snails had been collected during malacological surveys, individually checked for patent trematode infections by cercarial shedding and then preserved in 100% ethanol for molecular analysis [6]. Cercariae from infected *B. globosus* were preserved on Whatman FTA cards and identified using molecular methods as *S. haematobium* or *S. bovis* from two and five snails, respectively [6]. In addition, individual *B. globosus* and *B. nasutus* (also collected from Pemba), which were shedding two other trematode species, *Euclinostomum* sp. and *Echinostoma* sp. respectively (unpublished data), were tested to investigate assay specificity. Additionally, 94 *B. globosus* snails from Wambaa (Pemba) collected during November 2018, that were not shedding any trematode cercariae were tested for infections by PCR.

All samples that gave the 538 bp *Schistosoma* specific amplicon (Figure 7), were further subjected to the SIX PCR assay with the resulting amplicons purified and sequenced to confirm the species of the infection. *S. haematobium* and *S. bovis* species identity was confirmed by analysis of the four species SNP's that exist in the ITS2 region [7] between *S. haematobium* and *S. bovis* (Table 3).

#### 4.7. Testing the *Schistosoma cox1* rapid-diagnostic PCR (RD-PCR) for secondary species identification

The patent *B. globosus* snails collected from Pemba shedding either *S. haematobium* (n = 2) or *S. bovis* (n = 5) (see [6]), as detailed above, were further tested using the published multiplex RD-PCR (see [27,33]) with an aim to provide a secondary species-specific screening method as described in Schols *et al.* (2019) [27]. This multiplex RD-PCR, capable of differentiating *S. bovis* and *S. haematobium*

by species-specific amplicon size (*S. haematobium* (543 bp) *S. bovis* (306 bp), was performed following the published protocol and cycling conditions described by Webster *et al.* (2010) [33]. Different amount of gDNA (1 µl, 2 µl and 3 µl) and PCR annealing temperatures (58°C, 62°C and 65°C) were trialed to investigate sensitivity and specificity. Amplicons were purified and Sanger sequenced as described above, using the species-specific reverse primers to confirm species/amplicon identification.

#### 4.8. Ethical Statement

*Schistosoma haematobium* experimental infections were conducted at the Biomedical Research Institute – Schistosomiasis Resource Centre (Rockville, MA, USA) animal facility maintained with AAALAC full accreditation (Site # 000779), operating under the National Institutes of Health's Office of Laboratory Animal Welfare (OLAW) # A3080-01. *S. haematobium* parasite material was collected from male LVG Syrian golden hamsters following percutaneous exposure to cercariae. Hamster use was approved by the Institutional Animal Care and Use Committee (IACUC) of the Biomedical Research Institute for the Animal Use Protocol, #18-01.

**Author Contributions:** Conceptualization, Tom Pennance, Said M. Ali, Shaali M. Ame and Bonnie L. Webster; Data curation, Tom Pennance, Fiona Allan and Bonnie L. Webster; Formal analysis, Tom Pennance and Bonnie L. Webster; Funding acquisition, Bonnie L. Webster; Investigation, Tom Pennance, John Archer, Elena Lugli, Penny Rostron, Felix Llanwarne, Amour Khamis Amour, Khamis Rashid Suleiman, Sarah Li and Bonnie L. Webster; Methodology, Tom Pennance, John Archer, Sarah Li, Jo Cable and Bonnie L. Webster; Project administration, Tom Pennance, Said M. Ali, Shaali M. Ame and Bonnie L. Webster; Resources, Tom Pennance; Supervision, Said M. Ali, Jo Cable, Stefanie Knopp, Shaali M. Ame and Bonnie L. Webster; Validation, Tom Pennance, John Archer, David Rollinson, Fiona Allan and Bonnie L. Webster; Visualization, Tom Pennance; Writing – original draft, Tom Pennance; Writing – review & editing, John Archer, Elena Lugli, David Rollinson, Jo Cable, Stefanie Knopp, Fiona Allan and Bonnie L. Webster.

**Funding:** This work was funded by a Wellcome Trust Seed Award, grant number 207728. The APC was funded by Wellcome Trust Seed Award. T.P. was supported by the NERC GW4+ DTP and the Natural Environmental Research Council (NE/L002434/1). F.A. was supported by the Wellcome Trust on the SCAN project (104958/Z/14/Z).

**Conflicts of Interest:** The authors declare no conflict of interest. The funders had no role in the design of the study; in the collection, analyses, or interpretation of data; in the writing of the manuscript, or in the decision to publish the results.

#### References

- Colley, D.G.; Bustinduy, A.L.; Secor, W.E.; King, C.H. Human schistosomiasis. *Lancet* **2014**, *383*, 2253–2264.
- Brown, D.S. *Freshwater Snails of Africa and their Medical Importance*; 2nd ed.; Taylor & Francis: London, UK, 1994;
- WHO Schistosomiasis: Fact sheets Available online: <https://www.who.int/news-room/fact-sheets/detail/schistosomiasis> (accessed on Apr 27, 2020).
- De Bont, J.; Vercruysse, J. Schistosomiasis in cattle. *Adv. Parasitol.* **1998**, *41*, 285–364.
- De Bont, J.; Vercruysse, J. The epidemiology and control of cattle schistosomiasis. *Parasitol. Today* **1997**, *13*, 255–262.
- Pennance, T.; Ame, S.M.; Amour, A.K.; Suleiman, K.R.; Allan, F.; Rollinson, D.; Webster, B.L. Occurrence of *Schistosoma bovis* on Pemba Island, Zanzibar: implications for urogenital schistosomiasis transmission monitoring. *Parasitology* **2018**, *145*, 1727–1731.
- Pennance, T.; Allan, F.; Emery, A.; Rabone, M.; Cable, J.; Garba, A.D.; Hamidou, A.A.; Webster, J.P.; Rollinson, D.; Webster, B.L. Interactions between *Schistosoma haematobium* group species and their *Bulinus* spp. intermediate hosts along the Niger River Valley. *Parasit. Vectors* **2020**, *13*, 268.

8. Montresor, A.; Crompton, D.W.T.; Hall, A.; Bundy, D.A.P.; Savioli, L. *Guidelines for the evaluation of soil-transmitted helminthiasis and schistosomiasis at community level*; Geneva, 1998;
9. WHO *Ending the neglect to attain the Sustainable Development Goals – A road map for neglected tropical diseases 2021–2030*; Geneva, 2020;
10. Mutsaka-Makuvaza, M.J.; Matsena-Zingoni, Z.; Tshuma, C.; Ray, S.; Zhou, X.-N.; Webster, B.; Midzi, N. Reinfection of urogenital schistosomiasis in pre-school children in a highly endemic district in Northern Zimbabwe: a 12 months compliance study. *Infect. Dis. Poverty* **2018**, *7*, 102, doi:10.1186/s40249-018-0483-7.
11. Tchuente, L.-A.T.; Rollinson, D.; Stothard, J.R.; Molyneux, D. Moving from control to elimination of schistosomiasis in sub-Saharan Africa: time to change and adapt strategies. *Infect. Dis. Poverty* **2017**, *6*, 42.
12. King, C.H.; Sturrock, R.F.; Kariuki, H.C.; Hamburger, J. Transmission control for schistosomiasis—why it matters now. *Trends Parasitol.* **2006**, *22*, 575–582.
13. Rollinson, D.; Knopp, S.; Levitz, S.; Stothard, J.R.; Tchuente, L.A.; Garba, A.; Mohammed, K.A.; Schur, N.; Person, B.; Colley, D.G. Time to set the agenda for schistosomiasis elimination. *Acta Trop.* **2013**, *128*, 423–440.
14. Secor, W.; Colley, D. When Should the Emphasis on Schistosomiasis Control Move to Elimination? *Trop. Med. Infect. Dis.* **2018**, *3*, 85.
15. Pennance, T.; Person, B.; Muhsin, M.A.; Khamis, A.N.; Muhsin, J.; Khamis, I.S.; Mohammed, K.A.; Kabole, F.; Rollinson, D.; Knopp, S. Urogenital schistosomiasis transmission on Unguja Island, Zanzibar: characterisation of persistent hot-spots. *Parasit. Vectors* **2016**, *9*, 646.
16. Bayssade-Dufour, C. Chétotaxies cercariennes comparées de dix espèces de schistosomes. *Ann. Parasitol. Hum. Comparée* **1982**, *57*, 467–485.
17. Bayssade-Dufour, C.; Cabaret, J.; Ngendahayo, L.D.; Albaret, J.-L.; Carrat, C.; Chabaud, A.G. Identification of *Schistosoma haematobium*, *S. bovis* and *S. curassoni* by multivariate analysis of cercarial papillae indices. *Int. J. Parasitol.* **1989**, *19*, 839–846.
18. Cunningham, L.J.; Lingley, J.K.; Haines, L.R.; Ndung'u, J.M.; Torr, S.J.; Adams, E.R. Illuminating the Prevalence of *Trypanosoma brucei* s.l. in *Glossina* Using LAMP as a Tool for Xenomonitoring. *PLoS Negl. Trop. Dis.* **2016**, *10*, e0004441.
19. Garrod, G.; Adams, E.R.; Lingley, J.K.; Saldanha, I.; Torr, S.J.; Cunningham, L.J. Identification of *Trypanosoma brucei gambiense* and *T. b. rhodesiense* in vectors using multiplexed high-resolution melt analysis. *bioRxiv* **2020**, *In Review*, doi:10.1101/2020.04.21.052928.
20. Cook, D.A.N.; Pilotte, N.; Minetti, C.; Williams, S.A.; Reimer, L.J. A superhydrophobic cone to facilitate the xenomonitoring of filarial parasites, malaria, and trypanosomes using mosquito excreta/feces. *Gates Open Res.* **2017**, *1*, doi:10.12688/gatesopenres.12749.2.
21. Minetti, C.; LaCourse, J.E.; Reimer, L.; Stothard, J.R. Focusing nucleic acid-based molecular diagnostics and xenomonitoring approaches for human helminthiasis amenable to preventive chemotherapy. *Parasitol. Open* **2016**, *2*, e16.
22. Rathinasamy, V.; Hosking, C.; Tran, L.; Kelley, J.; Williamson, G.; Swan, J.; Elliott, T.; Rawlin, G.; Beddoe, T.; Spithill, T.W. Development of a multiplex quantitative PCR assay for detection and quantification of DNA from *Fasciola hepatica* and the intermediate snail host, *Austropeplea tomentosa*, in water samples. *Vet. Parasitol.* **2018**, *259*, 17–24.
23. Hamburger, J.; Hoffman, O.; Kariuki, H.C.; Muchiri, E.M.; Ouma, J.H.; Koech, D.K.; Sturrock, R.F.; King,



- C.H. Large-scale, polymerase chain reaction-based surveillance of *Schistosoma haematobium* DNA in snails from transmission sites in coastal Kenya: a new tool for studying the dynamics of snail infection. *Am. J. Trop. Med. Hyg.* **2004**, *71*, 765–773.
24. Lu, L.; Zhang, S.-M.; Mutuku, M.W.; Mkoji, G.M.; Loker, E.S. Relative compatibility of *Schistosoma mansoni* with *Biomphalaria sudanica* and *B. pfeifferi* from Kenya as assessed by PCR amplification of the *S. mansoni* ND5 gene in conjunction with traditional methods. *Parasites and Vectors* **2016**, *9*, 166.
25. Allan, F.; Rollinson, D.; Smith, J.E.; Dunn, A.M. Host choice and penetration by *Schistosoma haematobium* miracidia. *J. Helminthol.* **2009**, *83*, 33–38, doi:10.1017/S0022149X08073628.
26. Allan, F.; Dunn, A.M.; Emery, A.M.; Stothard, J.R.; Johnston, D.A.; Kane, R.A.; Khamis, A.N.; Mohammed, K.A.; Rollinson, D. Use of sentinel snails for the detection of *Schistosoma haematobium* transmission on Zanzibar and observations on transmission patterns. *Acta Trop.* **2013**, *128*, 234–240.
27. Schols, R.; Carolus, H.; Hammoud, C.; Mulero, S.; Mudavanhu, A.; Huyse, T. A rapid diagnostic multiplex PCR approach for xenomonitoring of human and animal schistosomiasis in a ‘One Health’ context. *Trans. R. Soc. Trop. Med. Hyg.* **2019**, *113*, 722–729.
28. Casotti, M.O.; Gryschek, R.C.B.; Paula, F.M. de; Gomes-Gouvêa, M.; Pinho, J.R.R.; Tuan, R.; Dias-Neto, E.; Luna, E.J. de A.; Espírito-Santo, M.C.C. do Molecular detection of prepatent *Schistosoma mansoni* infection in *Biomphalaria glabrata* snail vectors. *Rev. Inst. Med. Trop. Sao Paulo* **2020**, *62*.
29. Kaiglová, A.; Changoma, M.J.S.; Špajdelová, J.; Jakubcová, D.; Bírová, K. Urinary schistosomiasis in patients of rural medical health centers in Kwale county, Kenya. *Helminthologia* **2020**, *57*.
30. Abbasi, I.; Webster, B.L.; King, C.H.; Rollinson, D.; Hamburger, J. The substructure of three repetitive DNA regions of *Schistosoma haematobium* group species as a potential marker for species recognition and interbreeding detection. *Parasit. Vectors* **2017**, *10*, 364.
31. Kane, R.A.; Rollinson, D. Repetitive sequences in the ribosomal DNA internal transcribed spacer of *Schistosoma haematobium*, *Schistosoma intercalatum* and *Schistosoma matthei*. *Mol. Biochem. Parasitol.* **1994**, *63*, 153–156.
32. Webster, B.L. Isolation and preservation of schistosome eggs and larvae in RNeasy® facilitates genetic profiling of individuals. *Parasit. Vectors* **2009**, *2*, 50.
33. Webster, B.L.; Rollinson, D.; Stothard, J.R.; Huyse, T. Rapid diagnostic multiplex PCR (RD-PCR) to discriminate *Schistosoma haematobium* and *S. bovis*. *J. Helminthol.* **2010**, *84*, 107–114.
34. Schmaedick, M.A.; Koppel, A.L.; Pilotte, N.; Torres, M.; Williams, S.A.; Dobson, S.L.; Lammie, P.J.; Won, K.Y. Molecular xenomonitoring using mosquitoes to map lymphatic filariasis after mass drug administration in American Samoa. *PLoS Negl. Trop. Dis.* **2014**, *8*, e3087–e3087.
35. Pilotte, N.; Unnasch, T.R.; Williams, S.A. The Current Status of Molecular Xenomonitoring for Lymphatic Filariasis and Onchocerciasis. *Trends Parasitol.* **2017**, *33*, 788–798.
36. Pilotte, N.; Zaky, W.I.; Abrams, B.P.; Chadee, D.D.; Williams, S.A. A Novel Xenomonitoring Technique Using Mosquito Excreta/Feces for the Detection of Filarial Parasites and Malaria. *PLoS Negl. Trop. Dis.* **2016**, *10*, e0004641–e0004641.
37. Caron, Y.; Righi, S.; Lempereur, L.; Saegerman, C.; Losson, B. An optimized DNA extraction and multiplex PCR for the detection of *Fasciola* sp. in lymnaeid snails. *Vet. Parasitol.* **2011**, *178*, 93–99.
38. Cucher, M.A.; Carnevale, S.; Prepelitchi, L.; Labbé, J.H.; Wisnivesky-Colli, C. PCR diagnosis of *Fasciola hepatica* in field-collected *Lymnaea columella* and *Lymnaea viatrix* snails. *Vet. Parasitol.* **2006**, *137*, 74–82.
39. Kozak, M.; Wedrychowicz, H. The performance of a PCR assay for field studies on the prevalence of *Fasciola hepatica* infection in *Galba truncatula* intermediate host snails. *Vet. Parasitol.* **2010**, *168*, 25–30.

40. Magalhães, K.G.; Passos, L.K.J.; Carvalho, O. dos S. Detection of *Lymnaea columella* infection by *Fasciola hepatica* through Multiplex-PCR. *Mem. Inst. Oswaldo Cruz* **2004**, *99*, 421–424.
41. Mostafa, O.M.S.; Taha, H.A.; Ramadan, G. Diagnosis of *Fasciola gigantica* in snail using the polymerase chain reaction (PCR) assay. *Egypt. Soc. Parasitol.* **2003**, *33*, 733–742.
42. Velusamy, R.; Singh, B.P.; Raina, O.K. Detection of *Fasciola gigantica* infection in snails by polymerase chain reaction. *Vet. Parasitol.* **2004**, *120*, 85–90.
43. Kaplan, R.M.; Dame, J.B.; Reddy, G.R.; Courtney, C.H. A repetitive DNA probe for the sensitive detection of *Fasciola hepatica* infected snails. *Int. J. Parasitol.* **1995**, *25*, 601–610.
44. Carolus, H.; Muzarabani, K.C.; Hammoud, C.; Schols, R.; Volckaert, F.A.M.; Barson, M.; Huyse, T. A cascade of biological invasions and parasite spillback in man-made Lake Kariba. *Sci. Total Environ.* **2019**, *659*, 1283–1292.
45. Born-Torrijos, A.; Poulin, R.; Raga, J.A.; Holzer, A.S. Estimating trematode prevalence in snail hosts using a single-step duplex PCR: how badly does cercarial shedding underestimate infection rates? *Parasit. Vectors* **2014**, *7*, 243.
46. Kumagai, T.; Furushima-Shimogawara, R.; Ohmae, H.; Wang, T.-P.; Lu, S.; Chen, R.; Wen, L.; Ohta, N. Detection of Early and Single Infections of *Schistosoma japonicum* in the Intermediate Host Snail, *Oncomelania hupensis*, by PCR and Loop-Mediated Isothermal Amplification (LAMP) Assay. *Am. J. Trop. Med. Hyg.* **2010**, *83*, 542–548.
47. Tong, Q.-B.; Chen, R.; Zhang, Y.; Yang, G.-J.; Kumagai, T.; Furushima-Shimogawara, R.; Lou, D.; Yang, K.; Wen, L.-Y.; Lu, S.-H.; et al. A new surveillance and response tool: Risk map of infected *Oncomelania hupensis* detected by Loop-mediated isothermal amplification (LAMP) from pooled samples. *Acta Trop.* **2015**, *141*, 170–177.
48. Hamburger, J.; Ramzy, R.M.; Jourdane, J.; Ruppel, A. A polymerase chain reaction assay for detecting snails infected with bilharzia parasites (*Schistosoma mansoni*) from very early prepatency. *Am. J. Trop. Med. Hyg.* **1998**, *59*, 872–876.
49. Melo, F.L.; Gomes, A.L. do V.; Barbosa, C.S.; Werkhauser, R.P.; Abath, F.G.C. Development of molecular approaches for the identification of transmission sites of schistosomiasis. *Trans. R. Soc. Trop. Med. Hyg.* **2006**, *100*, 1049–1055.
50. Fernández-Soto, P.; Gandasegui Arahuetes, J.; Sánchez Hernández, A.; López Abán, J.; Vicente Santiago, B.; Muro, A. A loop-mediated isothermal amplification (LAMP) assay for early detection of *Schistosoma mansoni* in stool samples: a diagnostic approach in a murine model. *PLoS Negl. Trop. Dis.* **2014**, *8*, e3126–e3126.
51. Gandasegui, J.; Fernández-Soto, P.; Hernández-Goenaga, J.; López-Abán, J.; Vicente, B.; Muro, A. Biompha-LAMP: A New Rapid Loop-Mediated Isothermal Amplification Assay for Detecting *Schistosoma mansoni* in *Biomphalaria glabrata* Snail Host. *PLoS Negl. Trop. Dis.* **2016**, *10*, e0005225.
52. Abbasi, I.; King, C.H.; Muchiri, E.M.; Hamburger, J. Detection of *Schistosoma mansoni* and *Schistosoma haematobium* DNA by loop-mediated isothermal amplification: identification of infected snails from early prepatency. *Am. J. Trop. Med. Hyg.* **2010**, *83*, 427.
53. Caldeira, R.L.; Jannotti-Passos, L.K.; Dos Santos Carvalho, O. Use of Molecular Methods for the Rapid Mass Detection of *Schistosoma mansoni* (Platyhelminthes: Trematoda) in *Biomphalaria* spp. (Gastropoda: Planorbidae). *J. Trop. Med.* **2017**, *2017*, 8628971.
54. Kane, R.A.; Stothard, J.R.; Rollinson, D.; Leclipteux, T.; Evraerts, J.; Standley, C.J.; Allan, F.; Betson, M.; Kaba, R.; Mertens, P.; et al. Detection and quantification of schistosome DNA in freshwater snails using

- either fluorescent probes in real-time PCR or oligochromatographic dipstick assays targeting the ribosomal intergenic spacer. *Acta Trop.* **2013**, *128*, 241–249.
55. Hamburger, J.; He-Na; Abbasi, I.; Ramzy, R.M.; Jourdane, J.; Ruppel, A. Polymerase chain reaction assay based on a highly repeated sequence of *Schistosoma haematobium*: a potential tool for monitoring schistosome-infested water. *Am. J. Trop. Med. Hyg.* **2001**, *65*, 907–911.
  56. Rudko, S.P.; Reimink, R.L.; Froelich, K.; Gordy, M.A.; Blankespoor, C.L.; Hanington, P.C. Use of qPCR-Based Cercariometry to Assess Swimmer's Itch in Recreational Lakes. *Ecohealth* **2018**, *15*, 827–839.
  57. Rudko, S.P.; Reimink, R.L.; Peter, B.; White, J.; Hanington, P.C. Democratizing water monitoring: Implementation of a community-based qPCR monitoring program for recreational water hazards. *PLoS One* **2020**, *15*, e0229701.
  58. Rudko, S.P.; Turnbull, A.; Reimink, R.L.; Froelich, K.; Hanington, P.C. Species-specific qPCR assays allow for high-resolution population assessment of four species avian schistosome that cause swimmer's itch in recreational lakes. *Int. J. Parasitol. Parasites Wildl.* **2019**, *9*, 122–129.
  59. Webster, B.L.; Culverwell, C.L.; Khamis, I.S.; Mohammed, K.A.; Rollinson, D.; Stothard, J.R. DNA barcoding of *Schistosoma haematobium* on Zanzibar reveals substantial genetic diversity and two major phylogenetic groups. *Acta Trop.* **2013**, *128*, 206–217.
  60. Webster, B.L.; Emery, A.M.; Webster, J.P.; Gouvras, A.; Garba, A.; Diaw, O.; Seye, M.M.; Tchuente, L.A.T.; Simoonga, C.; Mwanga, J.; et al. Genetic diversity within *Schistosoma haematobium*: DNA barcoding reveals two distinct groups. *PLoS Negl. Trop. Dis.* **2012**, *6*, e1882.
  61. Emery, A.M.; Allan, F.E.; Rabone, M.E.; Rollinson, D. Schistosomiasis collection at NHM (SCAN). *Parasit. Vectors* **2012**, *5*, 185.
  62. Lewis, F.A.; Liang, Y.-S.; Raghavan, N.; Knight, M. The NIH-NIAID schistosomiasis resource center. *PLoS Negl. Trop. Dis.* **2008**, *2*, e267–e267, doi:10.1371/journal.pntd.0000267.

**Sample Availability:** Samples (DNA extracts of snails and parasites) are available from the authors upon appropriate request.



© 2020 by the authors. Submitted for possible open access publication under the terms and conditions of the Creative Commons Attribution (CC BY) license (<http://creativecommons.org/licenses/by/4.0/>).